



3D-Side

www.3Dside.eu

info@3Dside.eu

0032 10 81 35 48


CRANIOPLASTY

- Surgery date : August 2016
- Surgeon : Dr Vaz
- Hospital : Clinique universitaire St Luc
- **Challenge** : Frontal reconstruction of a skull defect with a **large surface area** (102 cm²)

Frontal Skull Defect of large surface area


P
L
A
N
N
I
N
G

The reconstruction was planned with our 3D design software



P
R
O
D
U
C
T
I
O
N

By means of high-tech 3D-Printing the **anatomical model** of the patient's skull and the **reconstruction mold** were produced



S
U
R
G
E
R
Y

During surgery we produced the **implant**, which perfectly fits the skull defect. This enabled recovering the patient's skull and restoring the face harmony

